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**DE**

**Operating manual**

**EN**

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**FR**

# ***PUREval MX-G***

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**1 General / Safety instructions**

1.1 General information .....	16
1.2 Field of use .....	16
1.3 Safety instructions .....	16
1.4 Safety regulations .....	17
1.5 Disclaimer of liability .....	17
1.6 Specific safety and work instructions .....	18

**2 System overview**

2.1 Structure of the ion exchange column .....	19
2.2 Technical data .....	20
2.2.1 Physical data .....	20
2.2.2 Versions .....	20

**3 Installation instructions / operation**

3.1 Setting up the glass column .....	21
3.2 Filling in the resin .....	21
3.3 Avoiding air inclusion / expelling air bubbles .....	21
3.4 Removal of the resins .....	21
3.5 Cleaning and care .....	22
3.6 Operation of the column according to the principle of gravity .....	22
3.7 Operation of the column with pump pressure .....	23

**4 Laboratory experiments using ion-exchange resins**

4.1 Experimental setup .....	24
4.2 Storage of the resins .....	24
4.3 Use of the resins for water-free applications .....	24
4.4 Rinsing .....	24
4.5 Pretreatment .....	24
4.6 Resin volume .....	25
4.7 Backwash .....	25
4.8 Recommended time and flow parameters .....	25
4.9 Sampling .....	26
4.10 Process parameters .....	26
4.11 Operation in laboratory tests .....	26
4.12 Regeneration .....	27
4.13 Advantages of the laboratory tests .....	27

## Chapter 1 - General / Safety instructions

### 1.1 General information

The ion-exchange column PURVEval MX-G is used to operate ion-exchange resins and adsorber resins on a laboratory scale.

The functional parts are made from high quality plastic parts. The ion exchange column consists of pressure-resistant borosilicate glass, all rubber parts are made of age-resistant elastomers. The materials used correspond to the recognized rules of technology.

Every person who works with this ion-exchange column must read these operating instructions in full beforehand and observe and apply the instructions given.

In addition to the operating instructions, the current and locally applicable regulations for accident prevention and for safe and professional work must be observed.

These operating instructions must always be available on site.

### 1.2 Field of use

The ion-exchange column is used to treat liquid process media using ion-exchange resins. The ion-exchange column is only suitable for technical purposes and not for the preparation of liquids for human consumption.

### 1.3 Safety instructions

Please read these operating instructions carefully before using the device and follow the instructions. The operating instructions must be kept ready to hand at all times.

Personal injury and property damage caused by non-observance of these operating instructions are not covered by the product liability law. The manufacturer assumes no liability for any other damage caused by non-compliance with these operating instructions.

Safety instructions warn of dangers and help to avoid personal injury and property damage. For your own safety, compliance with the safety instructions in these operating instructions is essential.

The applicable national and international safety regulations must be observed.

Each operator is responsible for compliance with the regulations that apply to him and must independently seek the latest regulations.

#### **1.4 Safety regulations**

The ion-exchange column may only be used by qualified personnel.

The specifications of the respective manufacturer must be observed for the maintenance or replacement of the consumables of the system.

The manufacturer's warranty expires if the device is modified.

The manufacturer accepts no liability for damage caused by improper commissioning.

The device must not be operated in potentially explosive areas. The device may only be put into operation if it is in perfect condition.

The device may only be used for the treatment of liquid process media. Check the ion-exchange column for possible damage before commissioning. The intended use within the performance limits must be ensured.

Before any repair work, the column must be disconnected from any media supply, any liquids must be properly disposed and the column must be rinsed with deionized water.

Damaged ion-exchange columns must be put out of operation immediately. Have defective or damaged columns repaired only by specialists authorized by the manufacturer. This is in your own interest. You thus prevent inadequate repairs and possible personal injury.

#### **1.5 Disclaimer of liability**

Use must be carried out exactly as described in this manual. The manufacturer is not liable for any damage, including consequential damage, that may arise from incorrect installation or incorrect use of the product.

### 1.6 Specific safety and work instructions

The PUREval MX-G ion-exchange column is only suitable for technical applications. Do not use the system to process liquids for human consumption.

Please observe the manufacturer-specific safety data sheets for the ion-exchange resin and for all acids and alkalis that you use to operate the system.

The borosilicate glass used for the ion exchange column corresponds to acid resistance class 1 in accordance with DIN 12 116 and alkali resistance class 2 in accordance with DIN ISO 695. Hydrofluoric acid, concentrated phosphoric acid and strong bases when temperatures > 100°C occur at the same time can damage the glass surface.

The ion-exchange column must not be opened or dismantled during operation. The column must not be opened without first releasing pressure.

Protect the column from direct sunlight and mechanical damage. Do not use near sources of heat or open flames.

When using drinking water from the pipe network, a safety fitting must be installed on the inlet connection in accordance with the European standard DIN EN 1717.

If the water pressure on the inlet side is greater than 2.5 bar, a pressure reducer must be installed in front of the ion exchange column.

To avoid contamination by particles, we recommend installing a filter <100 µm in front of the ion exchange column.

All parts must be installed in accordance with the country-specific guidelines. The ion exchange column must be operated within the specified ambient temperatures.

**Beware of frost:** after commissioning, storage and transport of water-filled systems below 4 ° C should be avoided. Frost can destroy the system.


The liquids treated with the ion-exchange column must be disposed of in accordance with the manufacturer's specifications and local regulations.

It is assumed that glass laboratory equipment is used in the normal manner in accordance with "good laboratory practice".

The operating personnel may only operate the ion-exchange column with suitable protective equipment.

## Chapter 2 - System overview

### 2.1 Structure of the ion exchange column

**2.2.2.1 Technical data**


**2.2.1 Physical data**

Operating pressure.....1,5 bar at 20°C  
 recommended operating / water temperature.....5 - 90°C  
 Ambient- and storage temperature  
 for column in use .....4 - 30°C  
 Water resistance class (DIN ISO 719).....Class 1  
 Acid class (DIN 12 116).....Class 1  
 Alkali class (DIN ISO 695).....Class 2  
 Hose connection – Ø outer .....10 mm  
 Drip nozzle – Ø outer .....10 mm  
 Flange opening – Ø inner .....25 mm  
 Operating direction .....vertical

Material.....borosilicate laboratory glass

**2.2.2 Versions**


Model	Rec. amount of resin	A (mm)	B (mm)	C (mm)	D (mm)	Articel no.
PUREval MX-G 100	100 ml	515	240	120	32	42154
PUREval MX-G 250	250 ml	615	340	170	43,4	42150
PUREval MX-G 500	500 ml	695	420	210	55,6	42151

## Chapter 3 - Installation instructions / operation

### 3.1 Setting up the glass column

The ion exchange column PUREval MX-G should only be operated in a vertical position. We recommend operating the glass column using a laboratory stand system with stand clamps.

### 3.2 Filling in the resins

Before the resins are filled in, the ion exchange column should already be attached vertically to a suitable bracket.

Close the regulating valve at the outlet of the column, dismantle the flange clamp and remove the upper part of the column.

Add water to the column and make sure that the resins filled in are completely covered with water. Using a suitable funnel, you can fill in the ion exchange resins through the flange opening with an inner diameter of 25 mm.

To avoid leaks, clean the gasket and sealing surface of the flange clamp.

### 3.3 Avoiding air inclusions / expelling air bubbles

To avoid air pockets in the resin bed, you should make sure that the resins are always covered with a sufficient amount of water. In addition, once filled resins should not run dry out, but should always be kept under liquid. To expel air bubbles, the column can be operated using the countercurrent principle with a low flow rate.

### 3.4 Removal of the resins

To remove the resins, you should first let the column run dry and rinse it with deionized water. Then dismantle the upper flange and remove the column from the stand clamp.

Mount a hose on the lower hose connection through which you can lead water. Rotate the column counter to normal operating direction so that the flange opening is facing down. Apply water pressure to the hose and flush the column with water. Catch the ion-exchange resins so that you can dispose it properly.

### 3.5 Cleaning and care

The filters and column should be cleaned of dust and dirt particles before first time use. To do this, hot hydrochloric acid and then distilled water are passed through the column and the filter plates.

The column should always be cleaned after use. If the filter plates are clean, cleaning the surfaces with distilled water is sufficient. If the filter plates are contaminated, chemical cleaning with solvents or a suitable rinsing solution should be carried out in the reverse operating direction. A flushing pressure of  $>2.5$  bar must be prevented. The column is then dried in the normal operating direction from top to bottom with clean air being blown through and additionally cleaned.

Rapidly changing temperatures must be avoided to avoid glass breakage.


Note: When using hot concentrated phosphoric acid, hydrofluoric acid and hot alkalis, an enlargement of the pores of the filter plates is unavoidable. These agents attack the glass surface and are therefore unsuitable as cleaning agents.

### 3.6 Operation of the column according to the principle of gravity

The laboratory tests with the ion-exchange column can be carried out without pressure according to the gravity principle.

The operation of the column is thus limited to direct current operation with direct current regeneration.

The following sketch shows an exemplary structure of column operation based on the principle of gravity.


### 3.7 Operation of the column with pump pressure

To simulate real process applications, the PUREval MX-G ion-exchange column allows a test setup with pump operation.

This means that different operating and regeneration variants can be simulated, such as direct current direct current, direct current counter current or floating bed operation.

The following sketch shows an example of the structure of column operation with pump operation.

## Chapter 4 - Laboratory experiments using ion-exchange resins

The laboratory test with ion-exchange resins is used to develop or optimize industrial processes. For this reason, the experimental operation of the laboratory ion-exchange system should simulate the real operating conditions as well as possible.

### 4.1 Experimental setup

Depending on the purpose of the experiment, the ion-exchange column can either be operated according to the gravity principle or with pumps.

If the ion-exchange column is operated according to the gravity principle, the supply hose should be laid in a U-shape so that the lowest point of the hose bend is just above the resin bed. This prevents the resin bed from drying out. The use of pumps allows a comfortable and realistic test procedure.

Regardless of the operating mode, sufficient process fluid, regeneration agent and fully demineralized water should be provided.

### 4.2 Storage of the resins

If the ion-exchange resins are to be stored before the experiment is carried out, seal the storage container tightly against the atmosphere and do not allow the resins to dry out. Avoid direct sunlight and warm or cold temperatures.

### 4.3 Use of the resins for water-free applications

If the resins are to be used in an anhydrous application area, the resins must be provided in dry form for best results. This avoids contamination of the treated product with water. Resin samples that you have received in hydrogenated form should be pretreated with suitable solvents such as acetone or alcohol, insofar as the specific process allows.

### 4.4 Rinsing

The resins should always be rinsed with deionized water before use in order to remove washable production and regeneration residues. In most applications, 5 - 10 bed volumes are sufficient.

### 4.5 Pretreatment

Before performing the laboratory test, ion-exchange resins should be pretreated so that complete hydrogenation and swelling of the polymer is ensured. Unless otherwise required for technical reasons, the polymer should always be completely swollen and soaked.

Ion-exchange resins should never be filled into a dry column. About half or a third of the column should be filled with deionized water before the pretreated ion-exchange resin is poured into the column.

Depending on the ion-exchange column used, the bed height should at least correspond to the position of the calibration mark.

#### 4.6 Resin volume

The PUREVAL MX-G ion-exchange columns are designed for different resin volumes depending on the version.

If the ionic load on the input liquid is very low, a large amount of test liquid is required until the resin bed is exhausted. This greatly increases the time required for each test run. In the case of low ionic loads, preliminary tests with smaller amounts of resin should be carried out if necessary.

#### 4.7 Backwash

As soon as the ion-exchange column is filled with resin, the resin should be backwashed using the countercurrent principle with deionized water for a period of 10-15 minutes. This process is known as classifying. The resin bed expands, larger particles arrange themselves below and light small particles rise upwards.

After backwashing, the resin bed should rest for about 5 minutes so that it can settle ideally.

After the backwashing process, the resin bed will increase somewhat depending on the particle size distribution. You should note the height of the resin bed or the bed volume (BV) for the following calculations.

In the next step the column should be emptied until there is only a 1 cm layer of water above the resin bed. Then you can begin to run the sample liquid into the column.

Use the regulating valve at the bottom of the column to regulate the flow rate through the resin bed.

#### 4.8 Recommended time and flow parameters

Parameter	
Resin volume	100 - 500 ml
Resin bed (minimal)	Depending on the type of column 120, 170 bzw. 210 mm.
Flow rate in normal operation	2 - 50 BV / h ( 8 - 20 BV / h typical)
Flow rate under regeneration	2 - 6 BV / h (2 BV / typical)
Regenerant contact time	15 - 60 minutes ( $\geq 30$ minutes are recommended)
Slow displacement flush	1 - 2 BV
Final rinse	2 - 10 BV

For the nature of the regeneration chemicals, their concentration and quantity, and other technical details, you should consult your resin manufacturer.

#### 4.9 Sampling

During the test, you should routinely test the resin in order to be able to make adjustments to the pH value if necessary. In treated water, the pH value can drop sharply, which has a strong influence on the alkalinity, breakthrough behavior, slip and resin capacity. Regular sampling of the treated liquid

also provides the necessary information about the loading of the column.

The samples should not be taken from the collecting container, but rather at the time of sampling from the drip nozzle of the ion exchange column. Following this principle you can record exact breakthrough curves.

#### 4.10 Process parameters

The ideal process conditions vary widely depending on the application. For processes that have not yet been fully developed or for which no empirical values exist, the optimal specific flow rate should first be determined as part of a laboratory test.

In most low load applications and conventional ion exchangers, the flow rate can be up to 50 BV / h. Even higher flow rates are sometimes used. In special process applications, particularly in applications with selective exchange resins, the flow rates should be at a level of 1 - 10 BV / h. In order to achieve maximum regeneration of the resins, the regeneration process is always carried out at low flow rates of 1 - 6 BV / h. The regeneration process is followed by a slow displacement rinse and a final rinse to ensure the best possible removal of the regenerant.

In water-free process applications, the resin is usually regenerated using water. In this case, the process medium must first be displaced with water. This process is called „sweetening off“ and after regeneration, when water is again displaced with the process medium, called „sweetening on“. These terms come from the sugar industry, where ion-exchange resins and absorber resins are used in a variety of ways.

#### 4.11 Operation in laboratory tests

Once the laboratory test has started, it should not be interrupted or paused, but should be carried out until the resins are completely exhausted. The experiment should not be paused. Most ion-exchange processes are reversible and once the flow is stopped the resin bed tries to reach a state of equilibrium. As soon as this happens, ions can move back from the resin into the liquid. This can cause the resins to run out prematurely and result in incorrect test results.

Under normal conditions the resin should always be covered with liquid. The column should never be emptied so far that air bubbles creep into the resin bed. Air bubbles are difficult to remove and worsen test results.

Three consecutive cycles should be performed with reproducible results before changing any operating conditions to optimize performance.

It typically takes two to three test cycles to get reliable results.

## 4.12 Regeneration

Many different processes have been developed for the regeneration of ion-exchange resins.

The simplest and most common method is DC regeneration. With this method, the regeneration agent and the process medium flow through the resin bed in the same direction.

The first step in cocurrent regeneration is always backwashing the resin with fully demineralized water in order to loosen the compaction of the resin bed and remove impurities. In general, the entry of particles into the resin bed should be prevented by prefiltration of the process solution.

After the backwash, the regeneration is carried out, followed by the slow displacement rinse and the final rinse.

A more effective regeneration of the resins can be achieved with countercurrent regeneration. With countercurrent regeneration, the regenerant flows through the resin bed against the normal operating direction, which results in significantly less slip. With the PUREval MX-G ion exchange column, countercurrent regeneration can be easily implemented with an additional pump.

In the case of countercurrent regeneration, which is normally carried out using the upflow principle, turbulence in the resins must be avoided. In order to stabilize the resin bed, the resins have to be driven to the upper filter plate at high speed. Alternatively, the free space above the resin bed can be filled with cotton wool or special inert resins.

In laboratory tests, deionized water is always used to dilute regeneration agents and for the rinsing processes. This simplifies the calculation of the capacity and also avoids additional loading of the resin bed with other ions.

## 4.13 Advantages of the laboratory tests

With controlled laboratory tests, a large amount of data and experience can be collected for individual process applications before an ion-exchange system is implemented in real operation.

The evaluation of all individual steps under controllable, reproducible conditions is the most reliable and economically sensible way to evaluate the behavior of the ion-exchange resins and to gather knowledge for a final ion-exchange system.

The laboratory tests are used to proactively optimize processes. With the data collected in the experiment, the concentration of regeneration agents, the flow rate and other decisive parameters can be optimized in advance before an ion-exchange system goes into operation on a real scale.

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elector GmbH